

Model PU-02

LITHIUM SHEET PUNCH

It is a handy tool for lithium sheet punching .
All parts are made of plastic material.

The product can be pushed out with a pusher
through the punch core if it sticks on the
punch end.

Rated thickness of the sheet is 0.5-3 mm; and
the rated diameter is 5-30mm. dimensions
other than standard spec can be ordered.
Please contact us to consult.



SPECIFICATION

- **Dimension:** 50W x 50D x 70H mm
- **Weight:** 300 Gram
- **Material:** POM (Overall)



Note: A test with clay sheet as the
bigger image. Upper image presents
the pusher.