

## Model PU-01

# Electrode/ Separator Punch

PU-01 is used preparing electrode and separator. It includes a mini press with collect tray and a punching mold.

The mold is mounted without any bolts. it means you can switch mold to different diameter very easily, just take-out and take-in .

The normal cutting diameter is from 6 to 24mm, and extra small or large size can be consulted.

Coin cell separator cutting should accompany a printer paper as explained in manual .



## SPECIFICATION

- **Dimensions**  
Press: 110W x 150D x 280H mm  
Mold: 40W x 50D x 70H mm
- **Weight**  
Whole: 4.3Kg  
Mold: 600g
- **Material:**  
Press Stand: Cast Iron  
Mold: Tool Steel

